**FOR INFORMATION ONLY.
WHEN PERFORMING
THE ASSAY ALWAYS REFER
TO PACKAGE INSERT
SUPPLIED
WITH THE KIT**

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For in vitro diagnostic use

GB

- For Professional Use Only
- Please refer to the U.S. Department of Health and Human Services (Bethesda, Md., USA) publication No. (CDC) 88-8395 on laboratory safety procedures or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Follow local guidelines for disposal of all waste material.

Caution

Material used in the preparation of human source reagent has been tested and found to be Non Reactive for HIV 1 and 2 Antibody, HCV Antibody and Hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they were potentially infectious.

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CanAg AFP EIA

Instructions for use

Enzyme immunometric assay kit
For 96 determinations

INTENDED USE

The CanAg AFP EIA kit is intended for the quantitative determination of α -Fetoprotein (AFP) in human serum.

SUMMARY AND EXPLANATION OF THE ASSAY

α -Fetoprotein (AFP), the foetal equivalent to albumin, is a 67 kDa glycoprotein produced during embryonic development and found in high concentrations in foetal serum and amniotic fluid. In normal non-pregnant adults AFP is present in low concentrations in serum. However AFP may be markedly increased in the serum from patients with cancer of the liver, testis or ovary. Quantitative determination of AFP in serum may be valuable in the management of patients with suspected or diagnosed liver cancer or germ cell tumours of the testis or ovary (1, 2).

PRINCIPLE OF THE TEST

The CanAg AFP EIA is a solid-phase, non-competitive immunoassay based upon the direct sandwich technique. Calibrators, controls and patient samples are incubated together with biotinylated Anti-AFP monoclonal antibody and horseradish peroxidase (HRP) labelled Anti-AFP monoclonal antibody in Streptavidin coated microstrips. After washing, buffered Substrate/ Chromogen reagent (hydrogen peroxide and 3, 3', 5, 5' tetra-methylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction a blue colour will develop if antigen is present. The intensity of the colour is proportional to the amount of AFP present in the samples.

The colour intensity is determined in a microplate spectrophotometer at 620 nm (or optionally at 405 nm after addition of Stop Solution). Calibration curves are constructed for each assay by plotting absorbance value versus the concentration for each calibrator. The AFP concentrations of patient samples are then read from the calibration curve.

REAGENTS

- Each CanAg AFP EIA kit contains reagents for 96 tests.
- The expiry date of the kit is stated on the label on the outside of the kit box.
- Do not use the kit beyond the expiry date.
- Do not mix reagents from different kit lots.
- Store the kit at 2–8°C. Do not freeze.
- Opened reagents are stable according to the table below provided they are not contaminated, stored in resealed original containers and handled as prescribed. Return to 2–8°C immediately after use.

| Component | Quantity | Storage and stability after first opening |
|-----------|----------|---|
|-----------|----------|---|

MICROPLA

| | | |
|-------------------|---------|---|
| Microplate | 1 Plate | 2–8°C until expiry date stated on the plate |
|-------------------|---------|---|

12 x 8 wells coated with Streptavidin. After opening, immediately return unused strips to the aluminium pouch, containing desiccant. Reseal carefully to keep dry.

| | | |
|------------------------|---------|---|
| AFP Calibrators | 6 vials | 2–8°C until expiry date stated on the vials |
|------------------------|---------|---|

| | | | | |
|-----|-----|-----|----------|-------------|
| CAL | AFP | 0 | 0 µg/L | 1 x 0.75 mL |
| CAL | AFP | 5 | 5 µg/L | 1 x 0.75 mL |
| CAL | AFP | 25 | 25 µg/L | 1 x 0.75 mL |
| CAL | AFP | 100 | 100 µg/L | 1 x 0.75 mL |
| CAL | AFP | 250 | 250 µg/L | 1 x 0.75 mL |
| CAL | AFP | 500 | 500 µg/L | 1 x 0.75 mL |

Human AFP in a Tris-HCl buffered salt solution containing bovine serum albumin, an inert yellow dye and 0.01% methyl-isothiazolone (MIT) as preservative. Ready for use.

| Component | Quantity | Storage and stability after first opening | | | |
|---|----------|---|---|-------------|--|
| AFP Controls | 2 vials | 2–8°C until expiry date stated on the vials | | | |
| <table border="1"><tr><td>CONTROL</td><td>AFP</td><td>1</td></tr></table> | CONTROL | AFP | 1 | 1 x 0.75 mL | |
| CONTROL | AFP | 1 | | | |
| <table border="1"><tr><td>CONTROL</td><td>AFP</td><td>2</td></tr></table> | CONTROL | AFP | 2 | 1 x 0.75 mL | |
| CONTROL | AFP | 2 | | | |

Human AFP in a Tris-HCl buffered salt solution containing bovine serum albumin, and 0.01% methyl-isothiazolone (MIT) as preservative. Ready for use.

| | | | | |
|---|-----------|--|--|--|
| <table border="1"><tr><td>BIOTIN</td><td>Anti-AFP</td></tr></table> | BIOTIN | Anti-AFP | | |
| BIOTIN | Anti-AFP | | | |
| Biotin Anti-AFP | 1 x 15 mL | 2–8°C until expiry date stated on the vial | | |

Biotin Anti-AFP monoclonal antibody from mouse, approximately 1 µg/mL. Contains phosphate buffered saline (pH 7.2), bovine serum albumin, bovine immunoglobulin, blocking agents, Tween 20, an inert blue dye and 0.01% methyl-isothiazolone (MIT) as preservative. To be mixed with Tracer, HRP Anti-AFP before use.

| | | | | |
|---|-------------|--|--|--|
| <table border="1"><tr><td>CONJ</td><td>Anti-AFP</td></tr></table> | CONJ | Anti-AFP | | |
| CONJ | Anti-AFP | | | |
| Tracer, HRP Anti-AFP | 1 x 0.75 mL | 2–8°C until expiry date stated on the vial | | |

Stock solution of HRP Anti-AFP monoclonal antibody from mouse, approximately 20 µg/mL. Contains preservatives. To be mixed with Biotin Anti-AFP before use.

| | | | | |
|--|-----------|--|--|--|
| <table border="1"><tr><td>SUBS</td><td>TMB</td></tr></table> | SUBS | TMB | | |
| SUBS | TMB | | | |
| TMB HRP-Substrate | 1 x 12 mL | 2–8°C until expiry date stated on the vial | | |

Contains buffered hydrogen peroxide and 3, 3', 5, 5' tetramethyl-benzidine (TMB). Ready for use.

| Component | Quantity | Storage and stability after first opening |
|-----------|----------|---|
|-----------|----------|---|

STOP

| | | |
|----------------------|-----------|--|
| STOP Solution | 1 x 15 mL | 2–8°C until expiry date stated on the vial |
|----------------------|-----------|--|

Contains 0.12 M hydrochloric acid. Ready for use.

WASHBUF 25X

| | | |
|-------------------------|-----------|--|
| Wash Concentrate | 1 x 50 mL | 2–8°C until expiry date stated on the bottle |
|-------------------------|-----------|--|

A Tris-HCl buffered salt solution with Tween 20. Contains Germall II as preservative. To be diluted with water 25 times before use.

Indications of instability

The TMB HRP-Substrate should be colourless or slightly bluish. A blue colour indicates that the reagent has been contaminated and should be discarded.

WARNINGS AND PRECAUTIONS

For in vitro diagnostic use.

- For professional use only.
- Please refer to the US Department of Health and Human Services (Bethesda, Md., US) publication No. (CDC) 88-8395 on laboratory safety or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Follow local guidelines for disposal of all waste material.

Caution

Material used in the preparation of human source reagent has been tested and found to be Non Reactive for HIV-1/2 Antibody, HCV Antibody and Hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they were potentially infectious.

SPECIMEN COLLECTION AND HANDLING

The CanAg AFP EIA is intended for use with serum. Collect blood by venipuncture and separate the serum according to common procedures. Samples can be stored at 2–8° C for 2 days. For longer periods it is recommended to store the samples at –20° C or below. Avoid repeated freezing and thawing of the samples. Allow frozen samples to thaw slowly, preferably at 2–8° C over night and then bring the samples to room temperature before analysis.

PROCEDURE

Materials required but not supplied with the kit

1. Microplate shaker

Shaking should be medium to vigorous. Longitudinal shaking approximately 200 strokes/min, oscillations 700-900/min.

2. Microplate wash device

Automatic plate wash capable of performing 1 and 6 washing cycles with a minimal fill volume of 350 µL/well/washcycle.

The Nunc Immuno-8 manual strip washer is recommended if an automatic microplatewash is not used.

3. Microplate spectrophotometer

With a wavelength of 620 nm and/or 405 nm and an absorbance range of 0 to 3.0.

4. Precision pipettes

With disposable plastic tips to deliver microlitre and millilitre volumes. An 8-channel pipette or respenser pipette with disposable plastic tips for delivery of 100 µL is useful but not essential.

5. Distilled or deionized water

For preparation of Wash Solution.

Procedural notes

1. A thorough understanding of this package insert is necessary to ensure proper use of the CanAg AFP EIA kit. The reagents supplied with the kit are intended for use as an integral unit. Do not mix identical reagents from kits having different lot numbers. Do not use the kit reagents after the expiry date printed on the outside of the kit box.
2. Reagents should be allowed to reach room temperature (20–25°C) prior to use. The assay should only be performed at temperatures between 20–25°C to obtain accurate results. Frozen specimens should be brought to room temperature slowly and must be gently but thoroughly mixed after thawing.
3. Before starting to pipette calibrators, controls and patient specimens it is advisable to mark the strips to be able to clearly identify the samples during and after the assay.
4. The requirement for efficient and thorough washing for separation of bound and unbound antigen and reagents from the solid-phase bound antibody-antigen complexes is one of the most important steps in an EIA. In order to ensure efficient washing make sure that all wells are completely filled to the top edge with wash solution during each wash cycle, that wash solution is dispensed at a good flow rate, that the aspiration of the wells between and after the wash cycles is complete and that the wells are empty. If there is liquid left, invert the plate and tap it carefully against absorbent paper.
 - Automatic strip washer: Follow the manufacturer's instructions for cleaning and maintenance diligently and wash the required number of wash cycles prior to and after each incubation step. It's highly recommended to use *strip* process mode and *overflow* wash mode with a dispensing volume of 800 µL. The aspiration/wash device should not be left standing with the Wash Solution for long periods, as the needles may get clogged resulting in poor liquid delivery and aspiration.
5. The TMB HRP-Substrate is very sensitive for contamination. For optimal stability of the TMB HRP-Substrate, pour the required amount from the vial to a carefully cleaned reservoir or preferably a disposable plastic tray to avoid contamination of the reagent. Be sure to use clean disposable plastic pipette tips (or respenser pipette tip).
6. Be sure to use clean disposable plastic pipette tips and a proper pipetting technique when handling samples and reagents. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid. A proper pipetting technique is of particular importance when handling the TMB HRP-Substrate Solution.

Protocol Sheet

CanAg AFP EIA REF **600-10**

Mix the components directly before use. Use shaking conditions according to the Instructions.

| Step | Bottle/Plate | Procedure | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|---------------------------|---|--|---------------|---------------------------------|----------------------|---|----|---|---|-----|---|---|-----|---|---|-----|---|---|-----|---|---|-----|---|---|-----|---|---|-----|---|---|-----|---|----|-----|----|
| 1. Prepare Wash Solution | WASHBUF 25X | Dilute 50 mL of Wash Concentrate with 1200 mL of distilled water or deionized water. | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | CONJ Anti-AFP | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Prepare Antibody Solution | BIOTIN Anti-AFP | Mix 50 μ L of Tracer, HRP Anti-AFP, with 1 mL of Biotin Anti-AFP per strip: | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | <table border="1"><thead><tr><th>No. of Strips</th><th>Tracer, HRP Anti-AFP (μL)</th><th>Biotin Anti-AFP (mL)</th></tr></thead><tbody><tr><td>1</td><td>50</td><td>1</td></tr><tr><td>2</td><td>100</td><td>2</td></tr><tr><td>3</td><td>150</td><td>3</td></tr><tr><td>4</td><td>200</td><td>4</td></tr><tr><td>5</td><td>250</td><td>5</td></tr><tr><td>6</td><td>300</td><td>6</td></tr><tr><td>7</td><td>350</td><td>7</td></tr><tr><td>8</td><td>400</td><td>8</td></tr><tr><td>9</td><td>450</td><td>9</td></tr><tr><td>10</td><td>500</td><td>10</td></tr></tbody></table> | No. of Strips | Tracer, HRP Anti-AFP (μ L) | Biotin Anti-AFP (mL) | 1 | 50 | 1 | 2 | 100 | 2 | 3 | 150 | 3 | 4 | 200 | 4 | 5 | 250 | 5 | 6 | 300 | 6 | 7 | 350 | 7 | 8 | 400 | 8 | 9 | 450 | 9 | 10 | 500 | 10 |
| No. of Strips | Tracer, HRP Anti-AFP (μ L) | Biotin Anti-AFP (mL) | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 1 | 50 | 1 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 2 | 100 | 2 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 3 | 150 | 3 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 4 | 200 | 4 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 5 | 250 | 5 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 6 | 300 | 6 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 7 | 350 | 7 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 8 | 400 | 8 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 9 | 450 | 9 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 10 | 500 | 10 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

| Preparation of reagents | Stability of prepared reagent |
|--|---|
| Wash Solution | 2 weeks at 2–25°C in a sealed container |
| Pour the 50 mL Wash Concentrate into a clean container and dilute 25- fold by adding 1200 mL of distilled or deionized water to give a buffered Wash Solution. | |

| | |
|---|------------------|
| Antibody Solution | 3 weeks at 2–8°C |
| Prepare the required quantity of Antibody Solution by mixing 50 µL of Tracer, HRP Anti-AFP with 1 mL of Biotin Anti-AFP per strip (see table below and the Protocol Sheet). | |

| No. of Strips | Tracer, HRP Anti-AFP (µL) | Biotin Anti-AFP (mL) |
|---------------|---------------------------|----------------------|
| 1 | 50 | 1 |
| 2 | 100 | 2 |
| 3 | 150 | 3 |
| 4 | 200 | 4 |
| 5 | 250 | 5 |
| 6 | 300 | 6 |
| 7 | 350 | 7 |
| 8 | 400 | 8 |
| 9 | 450 | 9 |
| 10 | 500 | 10 |
| 11 | 550 | 11 |
| 12 | 600 | 12 |

Be sure to use a clean plastic or glass bottle for preparation of the Antibody Solution.

Alternative: Pour the content of the Tracer, HRP Anti-AFP into the vial of Biotin Anti-AFP and mix gently. Make sure that all of the Tracer, HRP Anti-AFP is transferred to the vial of Biotin Anti-AFP.

NOTE: The Antibody Solution is stable for 3 weeks at 2–8°C. Do not prepare more Antibody Solution than will be used within this period and make sure that it is stored properly.

Assay procedure

Perform each determination in duplicate for calibrators, controls and patient samples. A calibration curve should be run with each assay. All reagents and samples must be brought to room temperature (20–25°C) before use.

1. Start to prepare Wash Solution and Antibody Solution. It is important to use clean containers. Follow the instructions carefully.
2. Transfer the required number of microplate strips to a strip frame. (Immediately return the remaining strips to the aluminium pouch containing a desiccant and reseal carefully). Wash each strip once with the Wash Solution. Do not wash more strips than can be handled within 30 min.
3. Pipette 25 µL of the AFP Calibrators (CAL 0, 5, 25, 100, 250, 500), controls (C1, C2) and patient samples (unknowns-Unk) into the strip wells according to the following scheme:

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 etc |
|---|------------|------------|-------|---|---|---|-------|
| A | Cal 0 | Cal 250 | Unk 1 | | | | |
| B | Cal 0 | Cal 250 | Unk 1 | | | | |
| C | Cal 5 | Cal 500 | Unk 2 | | | | |
| D | Cal 5 | Cal 500 | Unk 2 | | | | |
| E | Cal 25 | C1 | etc. | | | | |
| F | Cal 25 | C1 | | | | | |
| G | Cal 100 | C2 | | | | | |
| H | Cal 100 | C2 | | | | | |

4. Add 100 µL of Antibody Solution to each well using a 100 µL precision pipette (or an 8-channel 100 µL precision pipette). Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or the surface of the liquid.
5. Incubate the frame containing the strips for 1 hour (\pm 5 min) at room temperature (20–25°C) with constant shaking of the plate using a microplate shaker.
6. Wash each strip 6 times, using the wash procedure described in Procedural notes item 4.

7. Add 100 μ L of TMB HRP-Substrate to each well using the same pipetting procedure as in item 4. The TMB HRP-Substrate should be added to the wells as quickly as possible and the time between the addition to the first and last well should not exceed 5 min.
8. Incubate for 30 min (\pm 5 min) at room temperature with constant shaking. Avoid direct sunlight.
9. Immediately read the absorbance at 620 nm in a microplate spectrophotometer.

Option

If the laboratory does not have access to a microplate spectrophotometer capable of reading at 620 nm, the absorbance can be determined as follows:

- Alt. 9. Add 100 μ L of Stop Solution. Mix and read the absorbance at 405 nm in a microplate spectrophotometer within 15 minutes after addition of Stop Solution.

Measurement range

The CanAg AFP EIA measures concentrations between 0.5 and 500 μ g/L. If AFP concentrations above the measuring range are to be expected, it is recommended to dilute samples with normal human serum prior to analysis. **NOTE:** The serum used for dilution should also be measured in order to determine the endogenous AFP concentration (see "Calculation of results").

Quality control

AFP Control 1 and 2 may be used for validation of the assay series. Ranges of expected results are indicated on the vial labels. If values outside of the specified range are obtained, a complete check of reagents and reader performance should be made and the analysis repeated. Each laboratory may in addition prepare its own serum pools at different levels, which can be used as internal controls in order to assure the accuracy of the assay.

Reference material

The 1st International Standard IS 72/225 may be used as a reference standard. Values for AFP Calibrators and Controls were assigned against a set of in-house reference standards whose values are traceable to IS 72/225 using the conversion factor 0.83, i.e. 1 μ g/L corresponds to 0.83 kIU/L.

CALCULATION OF RESULTS

If a microplate spectrophotometer reader with built-in data calculation program is used, refer to the manual for the plate reader and create a program using the concentration stated on the labels of each of the AFP Calibrators.

For automatic calculation of AFP results it is recommended to use either of the following methods:

- Cubic spline curve fit method. Calibrator 0 should be included in the curve with the value 0 µg/L.
- Spline smoothed curve fit method. Calibrator 0 should be used as plate blank.
- Interpolation with point-to-point evaluation. Calibrator 0 should be included in the curve with the value 0 µg/L.
- Quadratic curve fit method. Calibrator 0 should be included in the curve with the value 0 µg/L.

NOTE: 4-parametric or linear regression should not be used.

For manual evaluation, a calibration curve is constructed by plotting the absorbance (A) values obtained for each AFP calibrator against the corresponding AFP concentration (in µg/L), see figure below. The unknown AFP concentrations can then be read from the calibration curve using the mean absorbance value of each patient specimen.

If samples in an initial analysis give AFP levels higher than 500 µg/L the samples should be diluted 1/10 and 1/100 with normal human serum and reanalyzed to obtain the accurate AFP concentration. **NOTE:** The sample used for dilution should also be measured in order to determine the endogenous AFP concentration.


The AFP concentration of the undiluted sample is then obtained as follows:

$$\text{Dilution 1/10: } 10 \times ([\text{AFP}]_{\text{Diluted sample}} - (0.9 \times [\text{AFP}]_{\text{Normal serum}}))$$

$$\text{Dilution 1/100: } 100 \times ([\text{AFP}]_{\text{Diluted sample}} - (0.99 \times [\text{AFP}]_{\text{Normal serum}}))$$

Example of results

| Specimen | | | Calibrator values | Mean abs value (A) | AFP (µg/L) |
|------------|-----|-----|-------------------|--------------------|------------|
| CAL | AFP | 0 | 0 µg/L | 0.036 | |
| CAL | AFP | 5 | 5 µg/L | 0.083 | |
| CAL | AFP | 25 | 25 µg/L | 0.282 | |
| CAL | AFP | 100 | 100 µg/L | 0.938 | |
| CAL | AFP | 250 | 250 µg/L | 1.914 | |
| CAL | AFP | 500 | 500 µg/L | 2.854 | |
| Specimen A | | | | 0.220 | 19.4 |
| Specimen B | | | | 1.686 | 208 |

Example (do not use this curve or table above to determine actual assay results).

LIMITATIONS OF THE PROCEDURE

The level of AFP cannot be used as absolute evidence for the presence or absence of malignant disease, and the AFP test should not be used in cancer screening. The results of the test should be interpreted only in conjunction with other investigations and procedures in the diagnosis of disease and the management of patients, and the AFP test should not replace any established clinical examination.

Anti-reagent antibodies (human anti-mouse antibody (HAMA) or heterophilic antibodies) in the patient sample may occasionally interfere with the assay, even though specific blocking agents are included in the buffer.

EXPECTED VALUES

AFP was measured in 40 apparently healthy men and 93 apparently healthy women. The mean value obtained was 2.8 µg/L with a standard deviation of 2.6. The lower and upper extremes of the normal range were examined using IFCC recommended non-parametric statistical treatment. The reference interval contains the central 95% fraction of the reference distribution. Reference limits may accordingly be estimated as the 2.5% (lower) and 97.5% (upper) fractiles. These limits cut off a fraction of 2.5% of the values in each tail of the reference distribution. Non-parametric estimates:

| Fraction | Reference limit (µg/L) | 90% confidence interval |
|----------------------------|------------------------|-------------------------|
| 2.5 th (lower) | 0.1 | 0.0–0.3 |
| 97.5 th (upper) | 10 | 8.7–14.6 |

It is recommended that each laboratory establish their own normal range to account for such local environmental factors as diet, climate, living conditions, patient selection, etc. As AFP levels have been shown to increase with age, the use of age-specific reference intervals have been suggested (2, 3). It should also be borne in mind that the individual patient's own baseline results provides the most important reference point for interpretation of marker results (3, 4).

Precision

Total precision was calculated according to NCCLS guideline EP5-A (6) using four levels of frozen pooled human serum containing added AFP and nine different CanAg AFP EIA reagent combinations. Each sample was randomly pipetted (n=2/analysis) and analysed twice each day over 20 days.

| Sample | Replicates | Mean (µg/L) | Within-run SD (µg/L) | Within-run CV % | Between-day SD (µg/L) | Between-day CV % |
|--------|------------|-------------|----------------------|-----------------|-----------------------|------------------|
| AFP 1 | 80 | 7.8 | 0.2 | 2.0 | 0.1 | 1.8 |
| AFP 2 | 80 | 23.2 | 0.4 | 1.8 | 0.3 | 1.4 |
| AFP 3 | 80 | 207 | 3.5 | 1.7 | 3.5 | 1.7 |
| AFP 4 | 80 | 416 | 6.6 | 1.6 | 8.5 | 2.0 |

Detection limit

The detection limit of the CanAg AFP EIA is ≤ 0.5 µg/L defined as the concentration corresponding to the mean of the absorbance values of the AFP calibrator 0 plus 2 standard deviations according to formula:

$$\frac{2 \times \text{SD CAL } 0}{\text{OD CAL } 5 - \text{OD CAL } 0} \times 5 \text{ } \mu\text{g/L}$$

Recovery

Spiked serum samples were prepared by adding aliquots of a sample with highly elevated AFP to normal serum samples. The recovery of the added antigen was in the range 90–106 %. **NOTE:** Recovery studies should **not** be performed using the kit calibrators.

Hook effect

No hook effect has been noticed for samples up to 40 000 µg/L. However, since patients with advanced hepatocellular carcinoma may show extremely high levels, false low results due to a high dose hook effect may be seen in specimens from these patients. In order to avoid reporting misleadingly low results due to a hook effect at higher concentrations, particularly in patients for whom markers are being measured for the first time, or when very high AFP values may be expected, it is recommended to assay specimens at two dilutions (i.e. neat and diluted 1:100 with normal human serum).

Linearity

Patient samples were diluted with normal human serum and analysed. The obtained values were within $\pm 10\%$ of the expected values.

Specificity

The CanAg AFP EIA is based on two mouse monoclonal antibodies, AFPK51 and AFPK57, targeting two separate antigenic determinants on the AFP molecule (5). The NCCLS guideline EP7-P (7) was followed to determine possible sources of interference. The following substances and concentrations were tested and found not to interfere with the test.

| | Concentration with no significant ($\pm 10\%$) interference |
|-------------------------|---|
| Lipemia (Intralipid®) | 10 mg/mL |
| Bilirubin, unconjugated | 0.6 mg/mL |
| Hemoglobin | 2 mg/mL |

WARRANTY

The performance data presented here were obtained using the assay procedure indicated. Any change or modification of the procedure not recommended by Fujirebio Diagnostics may affect the results, in which event Fujirebio Diagnostics disclaims all warranties expressed, implied or statutory including the implied warranty of merchantability and fitness for use.

REFERENCES

1. Johnson, P.J, (2002) Tumor Markers in Primary Malignancies of the liver. In "Tumor Markers: Physiology, pathobiology, technology and clinical applications", ed. Dimandis E.P. AACC Press, Washington pp 269-276.
2. Stenman, U-H and Alfthan, H. (2002) Markers for Testicular Cancer. In "Tumor Markers: Physiology, pathobiology, technology and clinical applications", ed. Dimandis E.P. AACC Press, Washington, pp 351-359.
3. Christiansen, M., Hogdall, C.K., Andersen, J.R. and Norgaard-Pedersen, B. (2001) Alpha-fetoprotein in plasma and serum of healthy adults: preanalytical, analytical and biological sources of variation and construction of age-dependent reference intervals. *Scand J Invest* 61: 205-216.
4. Trapé, J., Botargues, J.M., Porta, F., Ricós, C., Badal, J.M., Salinas, R., Sala, M., and Roca, A. (2003) Reference change value for α -Fetoprotein and its application in early detection of hepatocellular carcinoma in patients with hepatic disease. *Clin Chem* 49(7): 1209-1211.
5. Nustad, K., Paus, E., Kierulf, B., and Borner, O.P. (1998) Specificity and affinity of 30 monoclonal antibodies against Alpha-Fetoprotein. *Tumor Biol* 19: 293-300.
6. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
7. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).

CanAg® is a registered trademark of Fujirebio Diagnostics AB

Fujirebio Diagnostics AB
Elof Lindälvs gata 13
PO Box 121 32
SE-402 42 Göteborg
Sweden
Phone + 46 31 85 70 30
Fax + 46 31 85 70 40
info@fdab.com
www.fdab.com