


FUJIREBIO DIAGNOSTICS

Serum Products

Fujirebio Diagnostics, Inc. provides both specialty human and animal serum products used in manufacturing reagents, calibrators and controls. These high quality products are effective for use as calibrators and controls, and many can be tested with your own reagents to develop a custom pool that meets your objectives.

The following flowchart depicts the standard procedure used in the manufacture of our three basic types of specialty human serum products. Our processes are flexible such that each product can be modified to produce application specific product to meet your needs.



Human Serum Processing

Serum Products

Bulk Serum and Blood Products are manufactured with the same quality and manufacturing standard you trust.

- Bulk manufacturing
- Certified quality system (ISO 9001:2000 and ISO 13485:2003)
- Final 0.2 μm filtration to ensure a low bioburden product
- Custom pools can be manufactured to meet your specifications

Product Description	Catalog Number	Protein
BASE Human Serum	1016011	5.5 - 6.0 g/dL
BASE Human Serum	1016012	6.0 - 6.5 g/dL
CLARIFIED Human Serum	1016022	5.0 - 5.5 g/dL
CLARIFIED Human Serum	1016023	5.5 - 6.0 g/dL
CLARIFIED Human Serum	1016024	6.0 - 6.5 g/dL
CHARCOAL ADSORBED Human Serum	1016027	5.0 - 5.5 g/dL
CHARCOAL ADSORBED Human Serum	1016028	5.5 - 6.0 g/dL

American Biological Technologies, Inc.
940 Crossroads Blvd.
Sequin, TX 78155
USA

+1 800 531 7963 phone
+1 830 372 4130 fax
info@fdi.com e-mail

Fujirebio Diagnostics Inc.
201 Great Valley Parkway
Malvern, PA 19355-1307
USA

+1 610 240 3800 phone
+1 610 240 3803 fax
info@fdi.com e-mail
www.fdi.com website

